

Product datasheet for **AP02212SU-N**

Resistin (RETN) (64-88) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:2000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Resistin (aa 64-88) bTG- conjugated
Specificity:	This antibody detects Resistin peptide (aa 64-88); Resistin.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	resistin
Database Link:	Entrez Gene 56729 Human Q9HD89
Synonyms:	RETN, FIZZ3, HXCP1, RSTN, ADSF
Note:	LocusID 56729



[View online »](#)