

Product datasheet for **AP02213SU-S**

Secretin (SCT) (1-13) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:5.000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Secretin (aa 1-13) poly Lysin conjugated
Specificity:	This antibody detects Secretin (aa 1-13).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	secretin
Database Link:	Entrez Gene 6343 Human P09683
Synonyms:	secretin
Note:	Reference no.: LocusID 6343



[View online »](#)