

## Product datasheet for **AP02214SU-N**

### Secretin (SCT) (14-27) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:500).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Secretin (aa 14-27) poly Lysin conjugated
Specificity:	This antibody detects synthetic Secretin (aa 14-27).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	secretin
Database Link:	<a href="#">Entrez Gene 6343 Human P09683</a>
Synonyms:	secretin
Note:	LocusID 6343



[View online »](#)