

Product datasheet for **AP02215SU-S**

Secretory Immunglobuline A (sIgA) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:5000).
Reactivity:	Human
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Human sIgA
Specificity:	This antibody detects sIgA.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.



[View online »](#)