

Product datasheet for **AP02220SU-S**

DEFB4B (4-41) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:5000).
Reactivity:	Human
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Synthetic human β -Defensin 2 (aa 4-41)
Specificity:	This antibody detects synthetic β -Defensin 2 (aa 4-41).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	defensin beta 4B
Database Link:	Entrez Gene 100289462 Human O15263
Background:	This antibiotic peptide is locally regulated by inflammation. Defensins form a family of microbicidal and cytotoxic peptides made by neutrophils. Members of the defensin family are highly similar in protein sequence. hBD-2 is a antimicrobial peptide, which belongs to the distinct family of beta-defensins.
Synonyms:	Beta-defensin 4A, Beta-defensin 2, BD-2, hBD-2, DEFB102, DEFB2, Skin-antimicrobial peptide 1, DEFB4B, DEFB4A



[View online »](#)