

Product datasheet for **AP02224SU-S**

DEFB105A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:2.000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human β -Defensin 5 (aa 17-51)
Specificity:	This antibody detects beta-Defensin at aa 17-51.
Formulation:	State: Serum State: Liquid serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	defensin beta 105A
Database Link:	Entrez Gene 245908 Human Q8NG35
Background:	Beta 5 defensin has potential antibacterial activity, is specifically expressed in testis and belongs to the beta defensin family.
Synonyms:	Beta-defensin 105, DEFB105A, DEFB105B, DEFB5, BD5, Beta-defensin 5, DEFB-5, DEFB105



[View online »](#)