

Product datasheet for **AP02225SU-N**

DEFB106A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1/5.000)
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Beta-Defensin 6 (aa 1-38)
Specificity:	This antibody detects beta-Defensin.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	defensin beta 106A
Database Link:	Entrez Gene 245909 Human Q8N104
Synonyms:	Defensin beta 106, Beta-defensin 6, DEFB-6, DEFB6, DEFB106A, BD6, DEFB106B, DEFB106



[View online »](#)