

## Product datasheet for **AP02225SU-S**

### DEFB106A Rabbit Polyclonal Antibody

#### Product data:

|                        |   |
|------------------------|---|
| Product Type:          | Primary Antibodies  |
| Applications:          | ELISA   |
| Recommended Dilution:  | <b>ELISA</b> (1/5.000)  |
| Reactivity:            | Human   |
| Host:                  | Rabbit  |
| Clonality:             | Polyclonal  |
| Immunogen:             | Synthetic Human Beta-Defensin 6 (aa 1-38)   |
| Specificity:           | This antibody detects beta-Defensin.  |
| Formulation:           | State: Serum<br>State: Lyophilized Serum  |
| Reconstitution Method: | Restore in aqua bidest to initial volume.   |
| Conjugation:           | Unconjugated  |
| Storage:               | Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing. |
| Stability:             | Shelf life: One year from despatch.   |
| Gene Name:             | defensin beta 106A  |
| Database Link:         | <a href="#">Entrez Gene 245909 Human Q8N104</a>   |
| Synonyms:              | Defensin beta 106, Beta-defensin 6, DEFB-6, DEFB6, DEFB106A, BD6, DEFB106B, DEFB106         |



[View online »](#)