

## Product datasheet for **AP02228SU-N**

### Thymosin beta 10 (TMSB10) (1-12) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:1.000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Thymosin b10 (aa 1-12) BSA conjugated
Specificity:	This antibody detects Thymosin b10 (aa 1-12). There were no cross reactivities obtained with human Thymosin b4, Thymosin b9, Thymosin b15, Thymosin b4 peptide (aa 1-4), and Thymosin b4 (aa 1-14).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	thymosin beta 10
Database Link:	<a href="#">Entrez Gene 9168 Human P63313</a>
Synonyms:	TMSB10, PTMB10, THYB10, Thymosin beta-9



[View online »](#)