

Product datasheet for **AP02237SU-N**

Thymulin Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R
Recommended Dilution:	RIA. ELISA (1/15000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Thymulin
Specificity:	This antibody reacts with Thymulin. It does not react with Human Thymosin alpha 1 and Human Thymosin beta 4.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Background:	Thymulin, a nonapeptide secreted by the thymus, is an essential hormone for T lymphocyte differentiation and function.
Synonyms:	Thymic Factor



[View online »](#)