

Product datasheet for **AP02237SU-S**

Thymulin Rabbit Polyclonal Antibody

Product data:

| | |
|------------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, R |
| Recommended Dilution: | RIA. ELISA (1/15000). |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic Human Thymulin |
| Specificity: | This antibody reacts with Thymulin. It does not react with Human Thymosin alpha 1 and Human Thymosin beta 4. |
| Formulation: | State: Serum State: Lyophilized Serum |
| Reconstitution Method: | Restore in aqua bidest to initial volume. |
| Conjugation: | Unconjugated |
| Storage: | Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: One year from despatch. |
| Background: | Thymulin, a nonapeptide secreted by the thymus, is an essential hormone for T lymphocyte differentiation and function. |
| Synonyms: | Thymic Factor |



[View online »](#)