

## Product datasheet for **AP02239SU-S**

### Thyroglobulin (TG) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R
Recommended Dilution:	<b>IRMA.</b> <b>ELISA</b> (1/1000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Human Thyroglobulin
Specificity:	This antibody detects Thyroglobulin.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	thyroglobulin
Database Link:	<a href="#">Entrez Gene 7038 Human</a> <a href="#">P01266</a>
Synonyms:	Thyreoglobulin, TGN
Note:	LocusID 7038



[View online »](#)