

## Product datasheet for **AP02240PU-S**

### Thyroglobulin (TG) Mouse Monoclonal Antibody [Clone ID: ÜS3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	ÜS3
Applications:	ELISA, R
Recommended Dilution:	<b>IRMA.</b> <b>ELISA.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human Thyroglobulin
Specificity:	This antibody detects Thyroglobulin.
Formulation:	50 mM TRIS, pH 7.4 State: Purified State: Lyophilized purified IgG fraction Preservative: 0.02% Sodium Azide
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml.
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	thyroglobulin
Database Link:	<a href="#">Entrez Gene 7038 Human P01266</a>
Synonyms:	Thyreoglobulin, TGN



[View online »](#)