

Product datasheet for **AP02245SU-N**

Transglutaminase 2 (TGM2) Rabbit Polyclonal Antibody

Product data:

| | |
|------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | ELISA |
| Recommended Dilution: | ELISA (> 1:4.000). |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Recombinant human TissueTransglutaminase [TGc] |
| Specificity: | This antibody detects TissueTransglutaminase [TGc]. |
| Formulation: | State: Serum State: Lyophilized serum |
| Reconstitution Method: | Restore in aqua bidest to initial volume. |
| Conjugation: | Unconjugated |
| Storage: | Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: One year from despatch. |
| Gene Name: | transglutaminase 2 |
| Database Link: | Entrez Gene 7052 Human P21980 |
| Synonyms: | TGase C, TGase-H, Tissue transglutaminase |



[View online »](#)