

## Product datasheet for **AP02245SU-S**

### Transglutaminase 2 (TGM2) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (> 1:4.000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Recombinant human TissueTransglutaminase [TGc]
Specificity:	This antibody detects TissueTransglutaminase [TGc].
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	transglutaminase 2
Database Link:	<a href="#">Entrez Gene 7052 Human P21980</a>
Synonyms:	TGase C, TGase-H, Tissue transglutaminase



[View online »](#)