

Product datasheet for **AP02246SU-N**

TGF alpha (TGFA) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western blot.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Recombinant human Transforming Growth Factor a
Specificity:	This antibody detects Transforming Growth Factor a.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	transforming growth factor alpha
Database Link:	Entrez Gene 7039 Human P01135
Synonyms:	Transforming growth factor alpha, ETGF
Note:	LocusID 7039



[View online »](#)