

## Product datasheet for **AP02247SU-S**

### SP7 (220-234) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:100.000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Transcription Factor SP7 (aa 220-234)
Specificity:	This antibody detects Transcription Factor SP7 (aa 220-234).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	Sp7 transcription factor
Database Link:	<a href="#">Entrez Gene 121340 Human</a> <a href="#">Q8TDD2</a>
Synonyms:	Transcription factor Sp7, Osterix



[View online »](#)