

## Product datasheet for **AP02252SU-N**

### ANP (NPPA) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R
Recommended Dilution:	<b>RIA:</b> 1/3500 (See Ref.1-6). <b>ELISA:</b> 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Urodilatin (pro-ANP aa 95-101)- $\beta$ SAC KLH conjugated
Specificity:	This antibody detects Synthetic Human Urodilatin (aa 95-101), Human Urodilatin. There were no cross reactivities obtained with Human ANP (aa 1-28) (< 0.1%), Human pro-ANP (aa 1-98) (< 0.1%), Human pro-ANP (aa 79-98) (< 0.1%), Human BNP (< 0.1%) and Human CNP (< 0.1%).
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	natriuretic peptide A
Database Link:	<a href="#">Entrez Gene 4878 Human P01160</a>
Synonyms:	Atrial Natriuretic Peptide, Prepronatriodilatin, NPPA, ANP, PND
Note:	LocusID 4878



[View online »](#)