

## Product datasheet for **AP02256SU-S**

### Urokinase Plasminogen Activator-Receptor Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC
Recommended Dilution:	ELISA. Immunohistochemistry on frozen sections (1:50 - 1:100) and paraffin sections (1:50, recommended procedure for paraffin sections: pretreatment with 0.05 % Pronase E for 20 minutes).
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Purified human Urokinase Plasminogen Activator-Receptor
Specificity:	This antibody detects Urokinase Plasminogen Activator-Receptor. There were no cross reactivities obtained with human tPA and human PAI-2.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Note:	LocusID 5329



[View online »](#)