

Product datasheet for **AP02257SU-S**

Urotensin II (UTS2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:3000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Urotensin II, KLH conjugated
Specificity:	This antibody detects synthetic Urotensin II. There were no cross reactivities obtained with human Endothelin (aa 1-21) (< 0,1%); human big Endothelin (aa 1-38) (< 0,1%); human Urodilatin; human ANP (1-28) (< 0,1%); human BNP (aa 1-32) (< 0,1%); Adrenomedulin (aa 1-52) (< 0,1%); Relaxin (< 0,1%); and Somatostatin (< 0,1%).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	urotensin 2
Database Link:	Entrez Gene 10911 Human O95399
Synonyms:	Urotensin-II, U-II, UII, UNQ525/PRO1068
Note:	LocusID 10911



[View online »](#)