

## Product datasheet for **AP02259SU-N**

### VPAC2 (VIPR2) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human VIP-R2 (aa 3-14)
Specificity:	This antibody detects VIP-R2.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	vasoactive intestinal peptide receptor 2
Database Link:	<a href="#">Entrez Gene 7434 Human P41587</a>
Synonyms:	VIPR2, VIP2R, Vasoactive intestinal polypeptide receptor 2, PACAP-R-3, PACAPR3, VIPR-2, VIP-R2
Note:	LocusID 7434



[View online »](#)