

Product datasheet for **AP02265SU-S**

Vitamin D-binding protein Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Human GC Globulin
Specificity:	This antibody detects GC Globulin.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	GC, vitamin D binding protein
Database Link:	Entrez Gene 2638 Human P02774
Synonyms:	DBP, VDB, Group-specific component, Gc-globulin



[View online »](#)