

## Product datasheet for **AP02269SU-S**

### Aphrodisin (107-126) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western blot (1:1.000).
Reactivity:	Hamster
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Aphrodisin Peptide (aa107-126), polyLysine-conjugated
Specificity:	This antibody detects Aphrodisin.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Database Link:	<a href="#">Q9Z117</a>
Background:	Aphrodisin is a glycoprotein originally isolated from hamster vaginal discharge which was demonstrated to be involved in pheromonal effects on male hamsters. It selectively binds odorants and may be involved in carrying odorants within the nasal epithelium to putative olfactory receptors.



[View online »](#)