

Product datasheet for **AP02280SU-N**

RLN2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Porcine
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Porcine Relaxin
Specificity:	This antibody detects Relaxin.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Database Link:	<u>Entrez Gene 396891 Pig P01348</u>
Synonyms:	Relaxin A chain, Relaxin B chain



[View online »](#)