

## Product datasheet for **AP02281SU-N**

### RLN2 Guinea Pig Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Porcine
Host:	Guinea Pig
Clonality:	Polyclonal
Immunogen:	Porcine Relaxin
Specificity:	This antibody detects Relaxin.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Database Link:	<a href="#">Entrez Gene 396891 Pig P01348</a>
Synonyms:	Relaxin A chain, Relaxin B chain



[View online »](#)