

Product datasheet for **AP02288SU-N**

Gpnmb (551-568) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1/5000.
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Rat Osteoactivin peptide (aa551-568)
Specificity:	This antibody detects Osteoactivin peptide (aa551-568).
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	glycoprotein nmb
Database Link:	Entrez Gene 113955 Rat Q9QZF6
Synonyms:	Gpnmb



[View online »](#)