

Product datasheet for **AP02289SU-N**

Parathyroid Hormone (PTH) (1-34) Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	RIA (1:2.000).
Reactivity:	Rat
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Synthetic rat PTH aa (1-34) bTG- conjugated
Specificity:	This antibody detects PTH (aa 1-34). There were no cross reactivities obtained with PTH (aa 53-84).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	parathyroid hormone
Database Link:	P01270
Synonyms:	Parathormone, Parathyrin
Note:	LocusID 24694



[View online »](#)