

Product datasheet for **AP02290SU-N**

Parathyroid Hormone (PTH) (53-84) Rabbit Polyclonal Antibody

Product data:

| | |
|------------------------|------------------------------------------------------------------------------------------------|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, R |
| Recommended Dilution: | ELISA (1:3000). IRMA (1:500). |
| Reactivity: | Rat |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic rat PTH (aa 53-84) bTG-conjugated |
| Specificity: | This antibody detects PTH (aa 53-84). |
| Formulation: | State: Serum State: Lyophilized serum |
| Reconstitution Method: | Restore in aqua bidest to initial volume. |
| Conjugation: | Unconjugated |
| Storage: | Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: One year from despatch. |
| Gene Name: | parathyroid hormone |
| Database Link: | P01270 |
| Synonyms: | Parathormone, Parathyrin |



[View online »](#)