

## Product datasheet for **AP02290SU-S**

### Parathyroid Hormone (PTH) (53-84) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R
Recommended Dilution:	ELISA (1:3000). IRMA (1:500).
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic rat PTH (aa 53-84) bTG-conjugated
Specificity:	This antibody detects PTH (aa 53-84).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	parathyroid hormone
Database Link:	<a href="#">P01270</a>
Synonyms:	Parathormone, Parathyrin



[View online »](#)