

Product datasheet for **AP02291SU-N**

ANP (NPPA) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	RIA: 1/3,000 (See Ref.2-4).
Reactivity:	Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Rat Urodilatin / Rat ANP (aa 95-102) poly-Lysine conjugated
Specificity:	This antibody detects Rat and Mouse Urodilatin. There was no cross reactivity obtained with Rat BNP, Rat CNP, Rat ANP (aa 99-126), or Human Urodilatin.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	natriuretic peptide A
Database Link:	P01160
Synonyms:	Atrial Natriuretic Peptide, Prepronatriodilatin, NPPA, ANP, PND



[View online »](#)