

Product datasheet for **AP02298SU-N**

Prolactin (PRL) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, R, WB
Recommended Dilution:	RIA (1:160.000). Immunohistochemistry. Western blot.
Reactivity:	Sheep
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Sheep Prolactin
Specificity:	This antibody detects Prolactin.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	prolactin
Database Link:	P01236
Synonyms:	prolactin



[View online »](#)