

## Product datasheet for **AP02301SU-N**

### Vitamin D (Cholecalciferol) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Vitamin D
Specificity:	This antibody recognizes Vitamin D (Cholecalciferol).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore with aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Note:	CAS-No. 67-97-0



[View online »](#)