

## Product datasheet for **AP02329PU-N**

### JUND pSer255 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Suitable for use in Western blot (1:500~1:1000) and Immunohistochemistry on paraffin sections (1:50~1:100).
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized phosphopeptide derived from human JunD around the phosphorylation site of serine 255 (G-E-SP-P-P).
Specificity:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific phosphopeptide. The antibody against non-phosphopeptide was removed by chromatography using non-phosphopeptide corresponding to the phosphorylation site. JunD (phospho-Ser255) antibody detects endogenous levels of JunD only when phosphorylated at serine 255.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide and 50% glycerol State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	JunD proto-oncogene, AP-1 transcription factor subunit
Database Link:	<a href="#">Entrez Gene 3727 Human P17535</a>



[View online »](#)

**Background:**

JunD is the most broadly expressed member of the Jun family and the AP1 transcription factor complex. It has been found that primary fibroblasts lacking murine JunD displayed p53-dependent growth arrest, upregulated p19(ARF) expression, and premature senescence. In contrast, immortalized cell lines lacking JunD showed increased proliferation and higher cyclin D1 levels. These properties were reminiscent of the effects of oncogenic RAS expression on primary and established cell cultures. Furthermore, JunD *-/-* fibroblasts exhibited increased p53-dependent apoptosis upon ultraviolet irradiation and were sensitive to the cytotoxic effects of tumor necrosis factor- $\alpha$ . The antiapoptotic role of JunD was confirmed using an in vivo model of TNF-mediated hepatitis. The authors proposed that JunD protects cells from senescence, or apoptotic responses to stress stimuli, by acting as a modulator of the signaling pathways that link RAS to p53.

**Synonyms:**

JUND



**Product images:**

Figure 2. Western blot analysis of extracts from 293 cells using JunD (phospho-Ser255) antibody.



P-Peptide

-

+

Figure 1. Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using JunD (phospho-Ser255) antibody.