

# Product datasheet for AP02395PU-N

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Phospholipase C gamma 1 (PLCG1) pTyr783 Rabbit Polyclonal Antibody

**Product data:** 

**Product Type:** Primary Antibodies

Applications: WE

Recommended Dilution: Western blot: 1/500-1/1000.

Incubate membrane with diluted antibody in 5% nonfat milk, 1xTBS, 0.1% Tween-20 at 4°C

with gentle shaking, overnight.

Reactivity: Human, Mouse

**Host:** Rabbit

Clonality: Polyclonal

Immunogen: Peptide sequence around phosphorylation site of Tyrosine 783 (G-F-Ψ-V-E) derived from

Human PLC-gamma-1.

**Specificity:** This antibody detects endogenous levels of PLC-gamma 1 only when phosphorylated at

Tyrosine 783.

Formulation: PBS (without Mg2+ and Ca2+), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide and 50% Glycerol

State: Aff - Purified

State: Liquid purified IgG fraction

**Concentration:** lot specific

**Purification:** Affinity Chromatography using epitope-specific phosphopeptide. The antibody against

non-phosphopeptide was removed by chromatography using non-phosphopeptide

corresponding to the phosphorylation site

**Conjugation:** Unconjugated

**Storage:** Store the antibody (in aliquots) at -20°C.

Avoid repeated freezing and thawing.

**Stability:** Shelf life: One year from despatch.

Gene Name: phospholipase C gamma 1

Database Link: Entrez Gene 5335 Human

P19174





### Phospholipase C gamma 1 (PLCG1) pTyr783 Rabbit Polyclonal Antibody – AP02395PU-N

**Background:** Phospholipase C gamma (PLCG) binds to receptor tyrosine kinases, e.g. EGFR, PDGFR, leading

to phosphorylation of PLC-gamma at Tyrosine 771, 783 and 1245. Phosphorylation of Tyrosine 783 leads to enzymatic activation of PLC-gamma. PLC-gamma hydrolyzes phosphatidylinostitol-4,5-bisphosphate (PIP2), thus generating the second messengers

inostitol-1,4,5-triphosphate (IP3) and diacylglycerol (DAG).

Synonyms: PLC1, Phospholipase C-gamma-1, PLC-gamma-1, PLC gamma1 Phospholipase C-II, PLC-II, PLC-

148

Note: Molecular Weight: 155 kDa

### **Product images:**




Figure 1. Western blot analysis of extract from HL60 cells, untreated or treated with EGF using PLCG1 antibody (Phospho-Tyr783).