

Product datasheet for AP02447PU-S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PAK1 pThr423/402/421 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IF, IHC

Recommended Dilution: Immunofluorescence: 1/100-1/200.

Immunohistochemistry on Paraffin-Embedded Sections: 1/50-1/100.

Reactivity: Human, Mouse, Rat

Host: Rabbit

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized phosphopeptide derived from human PAK1

around the phosphorylation site of threonine 423 (R-S-TP-M-V).

Specificity: This antibody detects endogenous levels of PAK1/PAK2/PAK3 only when phosphorylated at

Threonine 423/402/421.

Formulation: PBS (without Mg2+ and Ca2+), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide and 50% Glycerol.

State: Aff - Purified

State: Liquid purified Ig fraction.

Concentration: lot specific

Purification: Affinity-chromatography using epitope-specific phosphopeptide. The antibody against non-

phosphopeptide was removed by chromatography using non-phosphopeptide corresponding

to the phosphorylation site.

Conjugation: Unconjugated

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: p21 (RAC1) activated kinase 1

Database Link: Entrez Gene 18479 MouseEntrez Gene 29431 RatEntrez Gene 5058 Human

Q13153

Synonyms: PAK 1, PAK-1, Alpha-PAK, PAK alpha, p21-activated kinase 1, p65-PAK, PAK 2, PAK-2, Gamma-

PAK, PAK gamma, PAK65, p21-activated kinase 2, p58, PAK 3, PAK-3, Beta-PAK, PAK beta, p21-

activated kinase 3, Oligophrenin-3, OPHN3





Protein Families: Druggable Genome, Protein Kinase, Stem cell - Pluripotency

Protein Pathways: Axon guidance, Chemokine signaling pathway, Epithelial cell signaling in Helicobacter pylori

infection, ErbB signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, MAPK signaling pathway, Natural killer cell mediated cytotoxicity, Regulation of actin cytoskeleton,

Renal cell carcinoma, T cell receptor signaling pathway

Product images:




Figure 1. Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue, using PAK1/PAK2/PAK3 (Phospho-Thr423/Thr402/Thr421) antibody.




Figure 2. Immunofluorescence staining of methanol-fixed HeLa cells using PAK1/PAK2/PAK3 (Phospho-Thr423/Thr402/Thr421) antibody (Red).