

## Product datasheet for **AP02460PU-N**

### H3FA (HIST1H3A) pSer10 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	<b>Western blot:</b> 1/500 - 1/1000; Incubate membrane with diluted antibody in 5% nonfat milk, 1X TBS, 0,1% Tween-20 at 4°C with gentle shaking, overnight. <b>Immunohistochemistry on paraffin-embedded sections:</b> 1/50 - 1/100. <b>Immunofluorescence:</b> 1/100 - 1/200.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic phosphopeptide derived from human Histone H3.1 around the phosphorylation site of serine 10 (R-K-S <sub>p</sub> -T-G).
Specificity:	Histone H3.1 (phospho-Ser10) antibody detects endogenous levels of Histone H3.1 only when phosphorylated at serine 10.
Formulation:	PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide and 50% Glycerol State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	histone cluster 1, H3a
Database Link:	<a href="#">P68431</a>



[View online »](#)

**Background:**

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fibre is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures.

Covalent modifications of the canonical core histones, including acetylation, phosphorylation, methylation, and monoubiquitination are used to mark nucleosomes to create chromatin domains with a range of functions. The information encoded by histone modifications can contribute to the formation and/or maintenance of transcriptionally active and inactive chromatin in response to various signalling pathways.

**Synonyms:**

H3/a, H3/b, H3/c, H3/d, H3/f, H3/h, H3/i, H3/j, H3/k, H3/l, HIST1H3A, H3FA, HIST1H3B, HIST1H3C, HIST1H3D, HIST1H3E, HIST1H3F, HIST1H3G, HIST1H3H, HIST1H3I, HIST1H3J


**Product images:**


Figure 2. Western blot analysis of extract from HeLa cells using Histone H3.1 antibody (Lane 1 and 2) and Histone H3.1 (phospho-Ser10) antibody AP02460PU (Lane 3 and 4).




Figure 4 Western Blot analysis of extracts from 293 cells untreated or treated with serum using Histone H3.1 (pSer10) antibody




Figure 1. Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Histone H3.1 (phospho-Ser10) antibody AP02460PU.

P-Peptide

-

+




Figure 3 Immunofluorescence staining of methanol-fixed HeLa cells using Histone H3.1 (pSer10) antibody