

## Product datasheet for **AP02555PU-S**

### **NF-κB p65 (RELA) Rabbit Polyclonal Antibody**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	IHC, WB
<b>Recommended Dilution:</b>	Western Blot: 1/500~1/1000. Immunohistochemistry: 1/50~1/100.
<b>Reactivity:</b>	Human, Mouse
<b>Host:</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	The antiserum was produced against synthesized non-phosphopeptide derived from human NF-κB p65 around the phosphorylation site of threonine 435 (E-G-TP-L-S).
<b>Specificity:</b>	NF-κB p65 antibody detects endogenous levels of total NF-κB p65 protein.
<b>Formulation:</b>	PBS(without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4 containing 150 mM NaCl, 0.02% Sodium Azide and 50% Glycerol. State: Aff - Purified State: Liquid purified IgG fraction.
<b>Concentration:</b>	lot specific
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store the antibody undiluted at -20°C. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>Gene Name:</b>	RELA proto-oncogene, NF-κB subunit
<b>Database Link:</b>	<a href="#">Entrez Gene 5970 Human Q04206</a>



[View online »](#)

**Background:**

NFκB (Nuclear Factor NF-kappa-B) is a pleiotropic transcription factor that plays a role in many biological processes, including inflammation, immunity, differentiation, cell growth, tumorigenesis, and apoptosis. It is found as a homo- or heterodimeric complex containing the Rel-like domain containing proteins NFκB p65 (RELA/p65), RELB, NFκB1/p105, NFκB1/p50, REL and NFκB2/p52. The heterodimeric NFκB p65/p50 complex is the most abundant one. The dimers bind to kappa-B sites at their target genes, with the affinity of the interaction dependent on the subunit composition of the dimer. Furthermore, different dimers act as transcriptional activators or repressors, with the NFκB p65/p50 and p65-c-Rel complexes acting as activators.


NFκB activity is controlled by several different mechanisms, including post-translational modifications, subcellular localisation and interactions with other coactivators or corepressors. NFκB complexes are held in the cytoplasm in an inactive state by interaction with members of the NFκB inhibitor (IκB) family. Typically, phosphorylation of IκB by IκB kinases (IKKs) in response to different activators leads to degradation of the inhibitor, allowing NFκB to translocate into the nucleus. The inhibitory effect of IκBs is primarily exerted through their interaction with NFκB p65.

NFκB p65 is ubiquitinated leading to its proteosomal degradation, which is required for termination of the NFκB response. Phosphorylation of NFκB p65 on S536 stimulates acetylation of K310 by CBP, enhancing transcriptional activity. NFκB p65 is also acetylated at K122, enhancing DNA binding and impairing the interaction with NFκBIA. The protein is deacetylated by HDAC3.


Invasion of a host by a pathogen is frequently associated with the activation of NF-κB, which coordinates various aspects of immune function required for resistance to infection.

**Synonyms:**

NF kappa B p65, NFκB p65, Transcription factor p65, Rel A, NFκB3

**Product images:**


Western blot analysis of extracts from COS7 cells using NF-κB p65 antibody (Line 1 and 2) and NF-κB p65 (pThr435) antibody (Line 3 and 4).



Peptide

-

+

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using NF-κB p65 antibody.