

Product datasheet for **AP03067PU-N**

Mouse IgG2b Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2b Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG2b paraproteins
Isotype:	IgG
Formulation:	Borate buffered saline pH8.2 with 0.09% Sodium Azide State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)