

Product datasheet for **AP05355HR-N**

Human IgG3 Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG3 Sheep Polyclonal Antibody
Applications:	ELISA, IHC
Recommended Dilution:	ELISA. Immunohistochemistry on Frozen Sections: 1/100-1/200.
Reactivity:	Human
Host:	Sheep
Immunogen:	Human IgG3 (kappa & lambda) purified from Human myeloma serum
Isotype:	IgG
Formulation:	PBS Label: HRP State: Liquid purified Ig fraction Preservative: 0.01% Thiomersal Label: Horseradish Peroxidase
Concentration:	0.5 mg/ml
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Note:	Antiserum Preparation Antiserum was raised by repeated immunisation of sheep with highly purified antigen.



[View online »](#)