

Product datasheet for AP05374FC-N

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

OriGene Technologies, Inc.

Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Hamster IgG (F(ab)2 specific) Goat Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Hamster IgG (F(ab)2 specific) Goat Polyclonal Antibody

Applications: FC

Recommended Dilution: Flow Cytometry: 1/100 - 1/200.

Reactivity: Hamster **Host:** Goat

Immunogen: Hamster IgG

Isotype: IgG

Formulation: Phosphate buffered saline containing 0.09% Sodium Azide, 0.5% Bovine Serum Albumin

Label: FITC

State: Liquid F(ab')2 fragmend

Label: Fluorescein Isothiocyanate Isomer 1

Concentration: lot specific

Purification: Affinity chromatography

Conjugation: FITC

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

This product is photosensitive and should be protected from light.