

Product datasheet for **AP05374PU-N**

Hamster IgG (F(ab)2 specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Hamster IgG (F(ab)2 specific) Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000-1/2000.
Reactivity:	Hamster, Reptiles
Host:	Goat
Immunogen:	Hamster IgG
Isotype:	IgG
Formulation:	PBS State: Purified State: Liquid purified F(ab')2 fragment Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography, Pepsin digestion
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)