

Product datasheet for **AP05375AC-N**

Mouse IgM Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgM Goat Polyclonal Antibody
Applications:	IF
Recommended Dilution:	Immunofluorescence.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgM paraproteins
Isotype:	IgG
Formulation:	Borate buffered saline, pH 8.2 containing 0.09% Sodium Azide (NaN ₃) , Stabilizing agent Label: APC State: Liquid purified IgM Label: Allophycocyanin
Purification:	Affinity chromatography
Conjugation:	APC
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light.



[View online »](#)