

## Product datasheet for **AP05376AP-N**

### Mouse IgE Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1:100 - 1:1000.
Reactivity:	Mouse
Host:	Goat
Isotype:	IgG
Formulation:	50mM HEPES pH7.1 containing 0.09% Sodium Azide and 0.2% Bovine Serum Albumin Label: AP State: Liquid purified IgG Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)