

## Product datasheet for **AP05376PU-N**

### Mouse IgE Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1:100 - 1:1000.
Reactivity:	Mouse
Host:	Goat
Isotype:	IgG
Formulation:	Phosphate buffered saline pH7.2 containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)