

Product datasheet for **AP05377PU-N**

Rat IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgA Goat Polyclonal Antibody
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA: 1µg/ml - 10µg/ml. Western Blot. Immunofluorescence.
Reactivity:	Reptiles, Rat
Host:	Goat
Isotype:	IgG
Formulation:	PBS pH7.2 containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)