

Product datasheet for **AP05378FC-N**

Rat IgG1 Goat Polyclonal Antibody

Product data:

| | |
|-----------------------|--|
| Product Type: | Secondary Antibodies |
| Product Name: | Rat IgG1 Goat Polyclonal Antibody |
| Applications: | FC |
| Recommended Dilution: | Flow Cytometry: Use 50 µl of 1/400 diluted antibody to label 1 x 10 ⁶ cells in 100 µl. |
| Reactivity: | Rat |
| Host: | Goat |
| Isotype: | IgG |
| Formulation: | Phosphate buffered saline pH 7.2 Label: FITC State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1 |
| Concentration: | lot specific |
| Purification: | Affinity Chromatography |
| Conjugation: | FITC |
| Storage: | Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. |



[View online »](#)