

Product datasheet for **AP05378HR-N**

Rat IgG1 Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG1 Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/10,000. Western Blot: 1/1,000-1/10,000.
Reactivity:	Rat
Host:	Goat
Isotype:	IgG
Formulation:	Phosphate buffered saline pH 7.2 Label: HRP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.1% Proclin-300 Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)