

## Product datasheet for AP05378PU-N

## **Rat IgG1 Goat Polyclonal Antibody**

**Product data:** 

**Product Type:** Secondary Antibodies

Product Name: Rat IgG1 Goat Polyclonal Antibody

**Applications:** ELISA, FC, WB

**Recommended Dilution: ELISA:** 1-10 μg/ml.

Western Blot. Flow Cytometry.

**Reactivity:** Rat

**Host:** Goat

**Immunogen:** Goats with highly purified antigen.

Isotype: IgG

**Formulation:** Phosphate buffered saline pH7.2 containing 0.09% Sodium Azide

State: Purified

State: Liquid purified IgG fraction

**Concentration:** lot specific

**Purification:** Affinity Chromatography

Conjugation: Unconjugated

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com