

Product datasheet for **AP05379FC-N**

Rat IgG2a Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2a Goat Polyclonal Antibody
Applications:	FC
Recommended Dilution:	Flow Cytometry: 1/500.
Reactivity:	Rat
Host:	Goat
Isotype:	IgG
Formulation:	Phosphate buffered saline pH7.2 containing 0.09% Sodium Azide and 0.2% Bovine Serum Albumin Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)