

## Product datasheet for AP05380PU-N

## Rat IgG2b Rabbit Polyclonal Antibody

**Product data:** 

**Product Type:** Secondary Antibodies

Product Name: Rat IgG2b Rabbit Polyclonal Antibody

**Applications:** ELISA, FC, WB

**Recommended Dilution:** ELISA: 1μg/ml - 10μg/ml.

Western Blot. Flow Cytometry.

Reactivity: Rat

**Host:** Rabbit

**Isotype:** lgG

**Formulation:** Phosphate buffered saline pH7.2 containing 0.09% Sodium Azide

State: Purified

State: Liquid purified IgG

**Concentration:** lot specific

**Purification:** Affinity chromatography

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com